



Product Datasheet

Histone H3 antibody (orb136531)


Description	Rabbit polyclonal antibody to Histone H3
Reactivity	All, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ChIP, IHC, IP, WB
Immunogen	A synthetic peptide of human HIST3H3(NP_003484.1).
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	17KDa
Uniprot ID	Q16695
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-200, IF/ICC: 1:50-200, IP: 1:50-200, ChIP: 1:50-200
Expiration Date	12 months from date of receipt.



IHC-P image of human lung cancer tissue ...



IHC-P staining of human liver cancer tis...



Immunohistochemical analysis of formalin...