



## Product Datasheet

### HIST2H2BE antibody (orb136534)


|                            |                                                              |
|----------------------------|--------------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to HIST2H2BE                      |
| <b>Species/Host</b>        | Rabbit                                                       |
| <b>Reactivity</b>          | Human, Mouse, Rat                                            |
| <b>Conjugation</b>         | Unconjugated                                                 |
| <b>Tested Applications</b> | IHC, WB                                                      |
| <b>Immunogen</b>           | Recombinant fusion protein of human HIST2H2BE (NP_003519.1). |
| <b>Preservatives</b>       | Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.   |
| <b>Concentration</b>       | 1 mg/ml                                                      |
| <b>Storage</b>             | Store at -20°C. Avoid freeze / thaw cycles.                  |
| <b>Note</b>                | For research use only                                        |
| <b>Isotype</b>             | IgG                                                          |
| <b>Clonality</b>           | Polyclonal                                                   |
| <b>Purity</b>              | Affinity purification                                        |
| <b>MW</b>                  | 14KDa                                                        |
| <b>Uniprot ID</b>          | <b>Q16778</b>                                                |
| <b>Entrez</b>              | <b>8349</b>                                                  |
| <b>Dilution Range</b>      | WB: 1:500-2000, IHC-P: 1:50-200                              |
| <b>Expiration Date</b>     | 12 months from date of receipt.                              |



Western blot analysis of extracts of var...



Immunohistochemical staining of human es...



Immunohistochemistry of paraffin-embedde...