



Product Datasheet

RAD50 antibody (orb136691)

Description	Rabbit polyclonal antibody to RAD50
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, IP, WB
Immunogen	A synthetic peptide of human Rad50(NP_005723.2).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	154KDa
Uniprot ID	Q92878
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-200, IF/ICC: 1:50-200
Expiration Date	12 months from date of receipt.



Immunoprecipitation analysis of 150ug ex...



Western blot analysis of extracts of var...