



Product Datasheet

IKZF1 antibody (orb136704)

Description	Rabbit polyclonal antibody to IKZF1
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ChIP, ICC, IF, WB
Immunogen	A synthetic peptide of human Ikaros(NP_006051.1).
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	50-70KDa
Uniprot ID	Q13422
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-200, IF/ICC: 1:50-200
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of BT4...



Immunofluorescence analysis of U2OS cell...