


## Product Datasheet

### COG5 Human Over-expression Lysate (orb1367314)

|                        |                                                                               |
|------------------------|-------------------------------------------------------------------------------|
| <b>Description</b>     | COG5 HEK293T cell transient overexpression lysate<br>(as WB positive control) |
| <b>Target</b>          | COG5                                                                          |
| <b>Note</b>            | For research use only                                                         |
| <b>Uniprot ID</b>      | <b>Q9UP83</b>                                                                 |
| <b>Expiration Date</b> | 6 months from date of receipt.                                                |



Western  
blot  
validation  
of Left: Cell  
ly...