



Product Datasheet

IGF2BP2 antibody (orb136741)

Description	Rabbit polyclonal to IGF2BP2
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, WB
Immunogen	A synthetic peptide of human IGF2BP2/IMP2(NP_006539.3).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	66kDa
Uniprot ID	Q9Y6M1
Dilution Range	WB: 1:200-2000, IHC-P: 1:50-100
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of Jur...



Immunohistochemical staining of mouse li...