


## Product Datasheet

### SUGT1 Human Over-expression Lysate (orb1367485)

|                        |                                                                             |
|------------------------|-----------------------------------------------------------------------------|
| <b>Description</b>     | SUGT1 HEK293T cell transient overexpression lysate (as WB positive control) |
| <b>Target</b>          | SUGT1                                                                       |
| <b>Note</b>            | For research use only                                                       |
| <b>Uniprot ID</b>      | <b>Q9Y2Z0, A8K7W3</b>                                                       |
| <b>Expiration Date</b> | 6 months from date of receipt.                                              |



Western blot validation of Left: Cell ly...