


Product Datasheet

ADAMTS5 antibody (orb136769)

Description	Rabbit polyclonal to ADAMTS5
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant fusion protein of human ADAMTS5 (NP_008969.2).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	73KDa
Uniprot ID	Q9UNA0
Entrez	11096
Dilution Range	WB 1:1000 - 1:5000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...