


Product Datasheet

CDC37 antibody (orb136773)

Description	Rabbit polyclonal to CDC37
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC, IP, WB
Immunogen	Recombinant fusion protein of human CDC37(NP_008996.1).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	50kDa
Uniprot ID	Q16543
Dilution Range	WB: 1:200-1000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...