



Product Datasheet

HP1 gamma antibody (orb136782)


Description	Rabbit polyclonal to HP1 gamma
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, IP, WB
Immunogen	Recombinant fusion protein of human HP1 gamma/HP1 gamma/CBX3(NP_009207.2).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	23kDa
Uniprot ID	Q13185
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human gastric tissue.



Western blot analysis of extract of various tissues.



Immunohistochemical staining of human lung tissue.