


Product Datasheet

EB1 antibody (orb136791)

Description	Rabbit polyclonal to EB1
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	Recombinant fusion protein of human MAPRE1(NP_036457.1).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	34kDa
Uniprot ID	Q15691
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-100, IP: 1:20-50
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of HeL...