


## Product Datasheet

### MON1B Human Over-expression Lysate (orb1368603)

|                        |                                                                             |
|------------------------|-----------------------------------------------------------------------------|
| <b>Description</b>     | MON1B HEK293T cell transient overexpression lysate (as WB positive control) |
| <b>Target</b>          | MON1B                                                                       |
| <b>Note</b>            | For research use only                                                       |
| <b>Uniprot ID</b>      | <b>Q6ZR87</b> , <b>Q7L1V2</b>                                               |
| <b>Expiration Date</b> | 6 months from date of receipt.                                              |



Western  
blot  
validation  
of Left: Cell  
ly...