



Product Datasheet

PDL1 antibody (orb136870)

Description	Rabbit polyclonal antibody to PDL1
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, WB
Immunogen	A synthetic peptide of human PD-L1/CD274(NP_054862.1).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	40-50KDa
Uniprot ID	Q9NZQ7
Dilution Range	WB: 1:500 - 2000, IF/ICC: 1:50 - 200
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...



Immunofluorescence analysis of RAW264.7 ...