


## Product Datasheet

### HS747E2A (C22orf31) Human Over-expression Lysate (orb1368792)

|                        |                                                                                |
|------------------------|--------------------------------------------------------------------------------|
| <b>Description</b>     | C22orf31 HEK293T cell transient overexpression lysate (as WB positive control) |
| <b>Target</b>          | C22orf31                                                                       |
| <b>Note</b>            | For research use only                                                          |
| <b>Uniprot ID</b>      | <b>O95567</b>                                                                  |
| <b>Expiration Date</b> | 6 months from date of receipt.                                                 |



Western blot validation of Left: Cell ly...