


## Product Datasheet

### WDR21A (DCAF4) Human Over-expression Lysate (orb1368885)

|                        |                                                                             |
|------------------------|-----------------------------------------------------------------------------|
| <b>Description</b>     | DCAF4 HEK293T cell transient overexpression lysate (as WB positive control) |
| <b>Target</b>          | DCAF4                                                                       |
| <b>Note</b>            | For research use only                                                       |
| <b>Uniprot ID</b>      | <b>B2RDD6, Q8WV16</b>                                                       |
| <b>Expiration Date</b> | 6 months from date of receipt.                                              |



Western blot validation of Left: Cell ly...