


## Product Datasheet

### RAB4B Human Over-expression Lysate (orb1369119)

|                        |                                                                             |
|------------------------|-----------------------------------------------------------------------------|
| <b>Description</b>     | RAB4B HEK293T cell transient overexpression lysate (as WB positive control) |
| <b>Target</b>          | RAB4B                                                                       |
| <b>Note</b>            | For research use only                                                       |
| <b>Uniprot ID</b>      | <b>P61018</b> , <b>A0A024R0K8</b>                                           |
| <b>Expiration Date</b> | 6 months from date of receipt.                                              |



Western  
blot  
validation  
of Left: Cell  
ly...