


Product Datasheet

CHRNA10 antibody (orb136953)

Description	Rabbit polyclonal to CHRNA10
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human CHRNA10(NP_065135.2).
Preservatives	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	45kDa
Uniprot ID	Q9GZZ6
Dilution Range	WB: 1:200-1000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...