


## Product Datasheet

# ZNF395 Human Over-expression Lysate (orb1369978)

|                        |                                                                              |
|------------------------|------------------------------------------------------------------------------|
| <b>Description</b>     | ZNF395 HEK293T cell transient overexpression lysate (as WB positive control) |
| <b>Target</b>          | ZNF395                                                                       |
| <b>Note</b>            | For research use only                                                        |
| <b>Uniprot ID</b>      | <b>Q7L9C8, Q9H8N7</b>                                                        |
| <b>Expiration Date</b> | 6 months from date of receipt.                                               |



Western  
blot  
validation  
of Left: Cell  
ly...