


Product Datasheet

FBXO32 antibody (orb137037)

Description	Rabbit polyclonal antibody to FBXO32
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human Fbx32/FBOX32 (NP_478136.1).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	40KDa
Uniprot ID	Q969P5
Entrez	114907
Dilution Range	WB: 1:1000-2000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...