


## Product Datasheet

### MTA3 Human Over-expression Lysate (orb1370394)

|                        |                                                                               |
|------------------------|-------------------------------------------------------------------------------|
| <b>Description</b>     | MTA3 HEK293T cell transient overexpression lysate<br>(as WB positive control) |
| <b>Target</b>          | MTA3                                                                          |
| <b>Note</b>            | For research use only                                                         |
| <b>Uniprot ID</b>      | <b>Q9BTC8</b>                                                                 |
| <b>Expiration Date</b> | 6 months from date of receipt.                                                |



Western  
blot  
validation  
of Left: Cell  
ly...