



Product Datasheet

LC3 antibody (orb13706)

Description	LC3 antibody
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	FC, IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human MAP1A heavy chain (2651-2750/3014 aa)
Target	MAP1A
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	326 kDa
Uniprot ID	Q9BXW4
Dilution Range	IHC-P=1:100-500, Flow-Cyt=1µg /test, IHC-F=1:100-500, IF=1:100-500
Expiration Date	12 months from date of receipt.



IHC-P of human glioma tissue (MAP1A anti...



Immunohistochemical analysis of formalin...