


Product Datasheet

CLOCK antibody (orb137114)

Description	Goat polyclonal to CLOCK
Species/Host	Goat
Reactivity	Bovine, Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	CLOCK (aa458-472)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Concentration	0.03-0.1 µg/ml.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	Approx 100kDa band observed in Mouse and Rat Skeletal Muscle lysates (calculated MW of 96.4kDa according to Mouse NP_031741.1 and 97.0kDa according to Rat NP_068628.1). Recommended concentration: 0.03-0.1 µg/ml.
Protein Sequence	DTSTPPRQHLPAAHEK
Clonality	Polyclonal
MW	95.3
Entrez	9575
Dilution Range	ELISA: 1:128000, WB: 0.03-0.1 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Mouse (A) and R...