



Product Datasheet

Asymmetric DiMethyl-Histone H3-R17 antibody (orb137154)


Description	Rabbit polyclonal antibody to Asymmetric DiMethyl-Histone H3-R17
Reactivity	All, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IP, WB
Immunogen	A synthetic asymmetric dimethylated peptide around R17 of human histone H3(NP_003520.1).
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
MW	17KDa
Uniprot ID	P68431
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-200, IF/ICC: 1:50-200, IP: 1:50-200
Expiration Date	12 months from date of receipt.



IHC-P image of human lung cancer tissue ...



IHC-P staining of human kidney cancer ti...



Immunohistochemical analysis of formalin...