



Product Datasheet

C7 antibody (orb137161)


Description	Goat polyclonal to C7
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	complement component 7
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 95kDa band observed in Human Liver lysates (calculated MW of 93.5kDa according to NP_000578.2). Recommended concentration: 0.1-0.3µg/ml. Primary incubation was 1 hour.. IHC: Paraffin embedded Human Placenta and Adrenal Gland. Recommended concentration: 5µg/ml.
Protein Sequence	NELTGQFRNRVINTK
Clonality	Polyclonal
MW	93.5
Entrez	730
Dilution Range	ELISA: 1:128000, WB: 0.1-0.3 µg/ml, IHC-P: 5 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human Liver lys...



Immunohistochemical staining of Human Pl...



Immunohistochemical staining of Human Ad...