



Product Datasheet

CREB3L2 antibody (orb137171)

Description	Goat polyclonal to CREB3L2
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	CREB3L2
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 26kDa band observed in Human Placenta lysates (calculated MW of 27.7kDa according to NP_001240704.1). Recommended concentration: 1-3µg/ml. IHC: In paraffin embedded Human Placenta shows staining in select nuclei. Recommended concentration, 5-10µg/ml.
Protein Sequence	DRKLSESEPGDGE
Clonality	Polyclonal
MW	57.4; 27.7
Entrez	64764
Dilution Range	ELISA: 1:128000, WB: 1-3 µg/ml, IHC-P: 5-10µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human Placenta ...



Immunohistochemical staining of Human Pl...