



Product Datasheet

HOXA5 antibody (orb137175)

Description	Goat polyclonal to HOXA5
Species/Host	Goat
Reactivity	Bovine, Human, Mouse, Porcine
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	HOXA5 (aa83-92)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 28kDa band observed in Human, Mouse, Rat and Pig Kidney lysates (calculated MW of 29.2kDa according to Mouse NP_034583.1). Recommended concentration: 1-3µg/ml.
Protein Sequence	EPRYSQPATS
Clonality	Polyclonal
MW	29.3
Entrez	3202
Dilution Range	ELISA: 1:8000, WB: 1-3 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of staining of Hum...



Western blot analysis of staining of Pig...