



Product Datasheet

Superoxide Dismutase antibody [Out of stock] (orb137215)

Description	Rabbit polyclonal to Superoxide Dismutase.
Reactivity	Bovine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	The superoxide dismutase enzyme from bovine erythrocytes (CAS # 9054-89-1).
Target	Superoxide Dismutase
Preservatives	PBS, 0.05% Sodium Azide, 0.05% rAlbumin
Concentration	0.5 mg/ml
Storage	Store at 4°C for short term applications. For long term storage, aliquot and store at -20°C.
Note	For research use only
Clonality	Polyclonal
Purity	Protein A purified
Dilution Range	IHC-P:10 µg/ml, WB:0.02-0.1 µg/ml
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of paraffin...



Western blot analysis of 1) human, 2) mo...