


## Product Datasheet

# LRP15 (LRRC3B) Human Over-expression Lysate (orb1372471)

|                        |                                                                              |
|------------------------|------------------------------------------------------------------------------|
| <b>Description</b>     | LRRC3B HEK293T cell transient overexpression lysate (as WB positive control) |
| <b>Target</b>          | LRRC3B                                                                       |
| <b>Note</b>            | For research use only                                                        |
| <b>Uniprot ID</b>      | <b>A1LNH5, Q96PB8</b>                                                        |
| <b>Expiration Date</b> | 6 months from date of receipt.                                               |



Western blot validation of Left: Cell ly...