


## Product Datasheet

### ITLN2 Human Over-expression Lysate (orb1372706)

|                        |                                                                                |
|------------------------|--------------------------------------------------------------------------------|
| <b>Description</b>     | ITLN2 HEK293T cell transient overexpression lysate<br>(as WB positive control) |
| <b>Target</b>          | ITLN2                                                                          |
| <b>Note</b>            | For research use only                                                          |
| <b>Uniprot ID</b>      | <b>Q8WWU7</b>                                                                  |
| <b>Expiration Date</b> | 6 months from date of receipt.                                                 |



Western  
blot  
validation  
of Left: Cell  
ly...