


## Product Datasheet

### BAX Human Over-expression Lysate (orb1373026)

|                        |                                                                           |
|------------------------|---------------------------------------------------------------------------|
| <b>Description</b>     | BAX HEK293T cell transient overexpression lysate (as WB positive control) |
| <b>Target</b>          | BAX                                                                       |
| <b>Note</b>            | For research use only                                                     |
| <b>Uniprot ID</b>      | <a href="#">Q07812</a>                                                    |
| <b>Expiration Date</b> | 6 months from date of receipt.                                            |



Western blot validation of Left: Cell ly...