


## Product Datasheet

### RG9MTD3 (TRMT10B) Human Over-expression Lysate (orb1373305)

|                        |                                                                               |
|------------------------|-------------------------------------------------------------------------------|
| <b>Description</b>     | TRMT10B HEK293T cell transient overexpression lysate (as WB positive control) |
| <b>Target</b>          | TRMT10B                                                                       |
| <b>Note</b>            | For research use only                                                         |
| <b>Uniprot ID</b>      | <b>Q6PF06</b>                                                                 |
| <b>Expiration Date</b> | 6 months from date of receipt.                                                |



Western blot validation of Left: Cell ly...