


Product Datasheet

NICE2 (S100A7A) Human Over-expression Lysate (orb1374643)

Description	S100A7A HEK293T cell transient overexpression lysate (as WB positive control)
Target	S100A7A
Form/Appearance	1 vial of 20 ug lyophilized gene specific transient over-expression cell lysate
Storage	The lysate can be shipped at ambient temperature. Upon receiving, store the sample at -20°C. Lysate samples can be reconstituted with SDS Sample Buffer. Avoid repeated freeze-thaw cycles after reconstitution. Lysate samples are stable for 12 months from date of receipt when stored at -20°C
Tag	C-Myc/DDK
Note	For research use only
Source	Expression Host: HEK293T
MW	Predicted Molecular Weight is 11.1 kDa
Uniprot ID	Q86SG5
Expiration Date	6 months from date of receipt.



Western blot validation of Left: Cell ly...