



Product Datasheet

VEGFR3 antibody (orb13753)


Description	VEGFR3 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	FC, IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human VEGFR-3 (901-1000/1298 aa)
Target	VEGFR3
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	151 kDa
Uniprot ID	P35916
Dilution Range	IHC-P=1:100-500, Flow-Cyt=1µg /test, WB=1:500-2000, IHC-F=1:100-500, ICC=1:100, IF=1:100-500, ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...



Flow cytometric analysis of A549 cells u...



Flow cytometric analysis of Splenocyte c...