


## Product Datasheet

### PADI6 Human Over-expression Lysate (orb1375943)

|                        |                                                                             |
|------------------------|-----------------------------------------------------------------------------|
| <b>Description</b>     | PADI6 HEK293T cell transient overexpression lysate (as WB positive control) |
| <b>Target</b>          | PADI6                                                                       |
| <b>Note</b>            | For research use only                                                       |
| <b>Uniprot ID</b>      | <b>Q6TGC4</b>                                                               |
| <b>Expiration Date</b> | 6 months from date of receipt.                                              |



Western  
blot  
analysis of  
Left: Cell  
lysa...