

Product Datasheet

CD22 antibody (orb137656)

Description

Mouse monoclonal to CD22. The HLA class I

whose products...

Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested

FC, ICC, IHC-Fr

Applications
Preservatives

Each vial contains 100 µl 1 mg/ml purified monoclonal antibody in PBS containing 0.09% sodium azide.

Form/Appearance

Each vial contains 100 µl 1 mg/ml purified monoclonal antibody in PBS containing 0.09% sodium azide.

Storage

Storage: The antibody is shipped at ambient temperature and may be stored at +4°C. For prolonged storage prepare appropriate aliquots and store at or below -20°C. Prior to use, an aliquot is thawed slowly in the dark at ambient temperature, spun down again and used to prepare working dilutions by adding sterile phosphate buffered saline (PBS, pH 7.2). Repeated thawing and freezing should be avoided. Working dilutions should be stored at +4°C, not refrozen, and preferably used the same day. If a slight precipitation occurs upon storage, this should be removed by centrifugation. It will not affect the performance or the concentration of the product

Note

For research use only

Isotype

IgG1

Clonality

Monoclonal

Source

3H4 is a mouse monoclonal IgG1 antibody derived by fusion of SP2/0-Ag14 mouse myeloma cells with spleen cells from BALB/c mice immunized with Daudi cells, a Burkitt's lymphoma cell line.

Uniprot ID

P35329

Expiration Date

12 months from date of receipt.