


## Product Datasheet

### ETL (ADGRL4) Human Over-expression Lysate (orb1376780)

|                        |                                                                              |
|------------------------|------------------------------------------------------------------------------|
| <b>Description</b>     | ADGRL4 HEK293T cell transient overexpression lysate (as WB positive control) |
| <b>Target</b>          | ADGRL4                                                                       |
| <b>Note</b>            | For research use only                                                        |
| <b>Uniprot ID</b>      | <b>Q9HBW9</b>                                                                |
| <b>Expiration Date</b> | 6 months from date of receipt.                                               |



Western blot analysis of Left: Cell lysa...