


## Product Datasheet

# SINHCAF Human Over-expression Lysate (orb1376835)

|                        |                                                                              |
|------------------------|------------------------------------------------------------------------------|
| <b>Description</b>     | FAM60A HEK293T cell transient overexpression lysate (as WB positive control) |
| <b>Target</b>          | SINHCAF                                                                      |
| <b>Note</b>            | For research use only                                                        |
| <b>Uniprot ID</b>      | <b>Q9NP50</b>                                                                |
| <b>Expiration Date</b> | 6 months from date of receipt.                                               |



Western  
blot  
analysis of  
Left: Cell  
lysate...