


## Product Datasheet

### MCAK (KIF2C) Human Over-expression Lysate (orb1377562)

|                        |                                                                             |
|------------------------|-----------------------------------------------------------------------------|
| <b>Description</b>     | KIF2C HEK293T cell transient overexpression lysate (as WB positive control) |
| <b>Target</b>          | KIF2C                                                                       |
| <b>Note</b>            | For research use only                                                       |
| <b>Uniprot ID</b>      | <b>A0A140VKF1, Q99661</b>                                                   |
| <b>Expiration Date</b> | 6 months from date of receipt.                                              |



Western blot analysis of Left: Cell lysa...