


## Product Datasheet

# HLADQA1 (HLA-DQA1) Human Over-expression Lysate (orb1378730)

|                        |                                                                           |
|------------------------|---------------------------------------------------------------------------|
| <b>Description</b>     | HLA HEK293T cell transient overexpression lysate (as WB positive control) |
| <b>Target</b>          | HLA-DQA1                                                                  |
| <b>Note</b>            | For research use only                                                     |
| <b>Uniprot ID</b>      | <b>Q8MH44, P01909, P01908, A0A173ADG5</b>                                 |
| <b>Expiration Date</b> | 6 months from date of receipt.                                            |



Western blot analysis of  
Left: Cell lysa...