

Product Datasheet

TRIF/TICAM1 Antibody (orb137880)

Description

TRIF/TICAM1 Antibody

Species/Host

Rabbit

Reactivity

Mouse

Conjugation

Unconjugated

Tested

WB

Applications

Immunogen

A synthetic peptide corresponding to a sequence at the C-terminus of mouse TRIF (693-708aa NNHMMWGHTGAQSDDK).

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

WB: The detection limit for TICAM1 is approximately 0.5ng/lane under reducing conditions. Tested Species: In-house tested species with positive results. Predicted Species: Species predicted to be fit for the product based on sequence similarities. Other applications have not been tested. Optimal dilutions should be determined by end users. . Add 0.2ml of distilled water will yield a concentration of 500ug/ml.

Isotype

Rabbit IgG

Clonality

Polyclonal

MW

79230 MW

Uniprot ID

[Q80UF7](#)

Dilution Range


Western blot, 0.1-0.5µg/ml, Mouse

Expiration Date

12 months from date of receipt.

1 2 3 4

250KD
130KD



Western blot
analysis of
Lane 1:
Recombi...