

Product Datasheet

Bonzo/CXCR6 Antibody (orb137892)

Description

Bonzo/CXCR6 Antibody

Species/Host

Rabbit

Reactivity

Mouse, Rat

Conjugation

Unconjugated

Tested

IHC, WB

Applications

Immunogen

A synthetic peptide corresponding to a sequence at the C-terminus of mouse Bonzo (331-351aa EDSSKTCSASHNVETTSMFQL), different from the related rat sequence by two amino acids, and from the related human sequence by three amino acids

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Tested Species: In-house tested species with positive results. Predicted Species: Species predicted to be fit for the product based on sequence similarities. By Heat: Boiling the paraffin sections in 10mM citrate buffer, pH6.0, for 20mins is required for the staining of formalin/paraffin sections. Other applications have not been tested. Optimal dilutions should be determined by end users. . Add 0.2ml of distilled water will yield a concentration of 500ug/ml.

Isotype

Rabbit IgG

Clonality

Polyclonal

MW

40468 MW

Uniprot ID


Q9EQ16

Dilution Range


Immunohistochemistry (Paraffin-embedded Section), 0.5-1µg/ml, Rat, Mouse, By HeatWestern blot, 0.1-0.5µg/ml, Mouse, Rat

Expiration Date

12 months from date of receipt.



Western blot analysis using Anti-CXCR6 a...



IHC(P) analysis of Rat Spleen Tissue usi...