


## Product Datasheet

### CUL4A Human Over-expression Lysate (orb1379080)

|                        |                                                                             |
|------------------------|-----------------------------------------------------------------------------|
| <b>Description</b>     | CUL4A HEK293T cell transient overexpression lysate (as WB positive control) |
| <b>Target</b>          | CUL4A                                                                       |
| <b>Note</b>            | For research use only                                                       |
| <b>Uniprot ID</b>      | <b>Q13619</b>                                                               |
| <b>Expiration Date</b> | 6 months from date of receipt.                                              |



Western  
blot  
analysis of  
Left: Cell  
lysa...