

Product Datasheet

Protein Wnt-3a Wnt3a Antibody (orb137923)

Description

Protein Wnt-3a Wnt3a Antibody

100KD —

Species/Host

Rabbit

55KD —

Reactivity

Human, Mouse

35KD —

Conjugation

Unconjugated

25KD —

Tested

WB

15KD —

Applications

Western blot
analysis
using Anti-
WNT3A a...

Immunogen

A synthetic peptide corresponding to a sequence at the C-terminus of human Wnt3a (252-267aa WVETLRPRYTIFYFKVPT), identical to the related mouse sequence.

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

WB: The detection limit for WNT3A is approximately 3ng/lane under reducing conditions. Tested Species: In-house tested species with positive results. Predicted Species: Species predicted to be fit for the product based on sequence similarities. Other applications have not been tested. Optimal dilutions should be determined by end users. . Add 0.2ml of distilled water will yield a concentration of 500ug/ml.

Isotype

Rabbit IgG

Clonality

Polyclonal

MW

39365 MW

Uniprot ID

P56704

Dilution Range

Western blot, 0.1-0.5µg/ml, Human, Mouse

Expiration Date

12 months from date of receipt.