

Product Datasheet

Angiotensinogen antibody (FITC) (orb13793)

Description	Angiotensinogen antibody (FITC)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	FITC
Tested Applications	IF
Immunogen	KLH conjugated synthetic peptide derived from human Angiotensin II type 1A receptor
Target	Angiotensin II type 1A receptor
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P30556
Expiration Date	12 months from date of receipt.