

Product Datasheet

EAAT3/SLC1A1 Antibody (orb137935)

Description

EAAT3/SLC1A1 Antibody

Species/Host

Rabbit

Reactivity

Mouse, Rat

Conjugation

Unconjugated

Tested

WB

Applications

Immunogen

A synthetic peptide corresponding to a sequence at the N-terminus of mouse EAAT3 (10-27aa DWRRFLRNHWLLSTVAA), identical to the related rat sequence.

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Tested Species: In-house tested species with positive results. Predicted Species: Species predicted to be fit for the product based on sequence similarities. Other applications have not been tested. Optimal dilutions should be determined by end users. . Add 0.2ml of distilled water will yield a concentration of 500ug/ml.

Isotype

Rabbit IgG

Clonality

Polyclonal

MW

56694 MW

Uniprot ID


P51906

Dilution Range

Western blot, 0.1-0.5µg/ml, Mouse, Rat

Expiration Date

12 months from date of receipt.



Western blot analysis of
Lane 1: Rat Liver
Liv...