

Product Datasheet

Monocarboxylate transporter 5 SLC16A4 Antibody (orb137961)

Description

Monocarboxylate transporter 5 SLC16A4



Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

IHC, WB

Applications

Immunohistochemical analysis of formalin...

Immunogen

A synthetic peptide corresponding to a sequence at the N-terminus of human SLC16A4 (134-148aa VVTTYFKKRLALST), different from the related rat and mouse sequences by two amino acids.

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage


Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.

Note


For research use only

Application notes

Tested Species: In-house tested species with positive results. Predicted Species: Species predicted to be fit for the product based on sequence similarities. By Heat: Boiling the paraffin sections in 10mM citrate buffer, pH6.0, for 20mins is required for the staining of formalin/paraffin sections. Other applications have not been tested. Optimal dilutions should be determined by end users. . Add 0.2ml of distilled water will yield a concentration of 500ug/ml.



Western blot analysis of Lane 1: Rat Tes...



Immunohistochemical analysis of formalin...

Isotype

Rabbit IgG

Clonality

Polyclonal

MW

54022 MW

Uniprot ID

015374

Dilution Range

Immunohistochemistry (Paraffin-embedded Section), 0.5-1µg/ml, Human, Rat, Mouse, By HeatWestern blot, 0.1-0.5µg/ml, Human, Rat, Mouse