


Product Datasheet

Rat CNTF ELISA Kit (orb138051)

Description	Rat CNTF ELISA Kit
Reactivity	Rat
Range	31.2 pg/ml - 2,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: A2-M200
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA, citrate).
MW	22854 MW
Uniprot ID	P20294
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is less than the lowest standard, 31.2 pg/ml.
Expiration Date	Please enquire.



Rat CNTF
PicoKine
ELISA Kit
standard
cur...