


Product Datasheet

Human TIMP-4 ELISA Kit (orb138055)

Description	Human TIMP-4 ELISA Kit
Reactivity	Human
Range	78 pg/ml - 5,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: C30-P224
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum, plasma (heparin, EDTA) and human milk.
MW	25503 MW
Uniprot ID	Q99727
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is around 1000 pg/ml.
Expiration Date	Please enquire.



Human
TIMP-4
PicoKine
ELISA Kit
standard...