


Product Datasheet

Mouse Prolactin ELISA Kit (orb138058)

Description	Mouse Prolactin ELISA Kit
Reactivity	Mouse
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: L32-C228
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, citrate).
MW	25496 MW
Uniprot ID	P06879
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:5, concentration in serum and plasma is 6-50 ng/ml.
Expiration Date	Please enquire.



Mouse
Prolactin
PicoKine
ELISA Kit
stand...