


Product Datasheet

Human TNFRSF17/BCMA ELISA Kit (orb138060)

Description	Human TNFRSF17/BCMA ELISA Kit
Reactivity	Human
Range	31.2 pg/ml - 2,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: M1-A54
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	20165 MW
Uniprot ID	Q02223
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:20, concentration in serum and plasma is 2.1-3 ng/ml.
Expiration Date	Please enquire.



Human
TNFRSF17/BCMA
PicoKine ELISA
Kit s...